





CÁMARAS DE CLIMA CONSTANTE KBF

Modelo	Volumen Efectivo (L)	Dimensiones externas (bxaxp) cm	Dimensiones internas (bxaxp) cm	Rango de Temperatura (°C)	Intervalo de Humedad (H.R.)
 KBF 1020	1020	125.0x192.5x114.5	97.3x125.0x83.6	0 a 70	10% a 80%
 KBF 720	700	125.0x192.5x89.0	97.3x125.0x57.6	0 a 70	10% a 80%
 KBF 240	247	92.5x146.0x80.0	65.0x78.5x48.5	0 a 70	10% a 80%
 KBF 115	102	88.0x105.0x65.0	60.0x48.3x35.1	0 a 70	10% a 80%