

CÁMARAS DE CLIMA VARIABLE MKF

Modelo	Volumen Efectivo (L)	Dimensiones externas (bxaxp) cm	Dimensiones internas (bxaxp) cm	Rango de Temperatura (°C)	Intervalo de Humedad (H.R.)
 MKF 1020	991	158.0x200.5x148.0	120.0x102.0x81.0	-40 a 180	10% a 98%
 MKF 720	734	158.0x200.5x114.0	120.0x102.0x60.0	-40 a 180	10% a 98%
 MKF 400	417	111.5x171.0x140.0	73.5x70.0x81.0	-40 a 180	10% a 98%
 MKF 240	228	111.5x171.5x92.5	73.5x70.0x44.3	-40 a 180	10% a 98%
 MKF 115	115	98.0x172.5x86.5	60.0x48.0x40.0	-40 a 180	10% a 98%
 MKF 56	60	72.0x144.5x78.0	40.0x42.0x34.8	-40 a 180	10% a 98%