

CÁMARA DE CALOR Y SECADO ED

Modelo	Volumen Efectivo (L)	Dimensiones externas (bxaxp) cm	Dimensiones internas (bxaxp) cm	Rango de Temperatura (°C)	
	FD 720	741	116.5x159.0x87.0	100.0x130.0x57.0	5 a 300 Sobre la temperatura ambiente
	FD 260	259	81.0x94.0x82.5	65.0x78.0x51.5	5 a 300 Sobre la temperatura ambiente
	FD 115	116	71.0x73.5x60.5	55.0x55.0x38.5	5 a 300 Sobre la temperatura ambiente
	FD 56	60	56.0x62.5x56.5	40.0x44.0x34.5	5 a 300 Sobre la temperatura ambiente
	FD 23	20	43.5x49.5x52.0	22.2x33.0x30.0	5 a 300 Sobre la temperatura ambiente