






INCUBADORA REFRIGERADA KB

Modelo	Volumen Efectivo (L)	Dimensiones externas (bxaxp) cm	Dimensiones internas (bxaxp) cm	Rango de Temperatura (°C)	
	KB 720	698	125.0x192.5x88.5	97.0x125.0x57.6	-10 a 110
	KB 400	400	92.5x195.0x80.5	65.0x127.0x48.5	-10 a 110
	KB 240	247	92.5x146.5x80.0	65.0x78.5x48.5	-10 a 110
	KB 115	115	83.5x102.5x65.0	60.0x48.0x40.0	-10 a 110
	KB 53	53	63.5x83.5x58.0	40.0x40.0x33.0	-10 a 110